

ROTARY EVAPORATOR OVERVIEW



RE200/210

- Manual lift with analog control
- Option for 3 different types of glassware
- Option for either a 4L or 7L capacity water bath
- Option for 7L capacity oil bath
- Optional arm jack permits setting flexibility



RE301

- Easy manual operation
- Digitally set and displayed motor speed (rpm)
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC



RE601

- Programmable
- Digital setting and display for motor speed (rpm), vacuum control and vapor temperature
- Equipped with vacuum regulator and vapor temperature indicator
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC



RE801

- Programmable
- Bath synchronized control
- Automatic distillation and 53 solvent presets
- Digital setting and display for motor speed (rpm), vacuum control and vapor temperature
- Equipped with vacuum regulator and vapor temperature indicator
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC