


## Stands


### Stands and suitable oven models

Product code	Stand models	Suitable oven models
211856	ON61	DX 400/410/401/402/402C/412C/600/601/602/610/602C/612C, DVS 400/401/402/402C/412C/600/601/602/611/612/602C/612C, DG400 DKM 300C/310C/400C/410C/600C/610C, DKN 300/301/302/302C/312C/400/401/402/402C/412C/600/601/602/611/612/602C/612C, DNE 400/401/410/411/600/610/611, DNF 400/410/600/610
212802	ONS60	DR200
415464	OP43 (without caster)	DF 411/412, DH 411/412
415465	OP63 (without caster)	DF 611/612, DH 611/612
415466	OP46 (with caster, front side with two stoppers)	DF 411/412, DH 411/412
415467	OP66 (with caster, front side with two stoppers)	DF 611/612, DH 611/612


#### ON61



#### OP43/63, OP46/66 (casters front side with two stoppers)




#### ONS60



## Shelves

### Stainless steel shelf & bracket set and suitable oven models


Product code	Punching shape	Suitable oven models
212246	Round punch shelf & bracket set	DVS 402/402C/412C, DG400, DKN 402/402C/412C, DNE 400/401/410/411, DNF 400/410
212266	Round punch shelf & bracket set	DX 600/601/602/610/602C/612C, DVS 600/601/602/611/612/602C/612C, DKN 600/601/602/611/612/602C/612C/810/811/812/812C, DNE 600/610/611/810/811, DNF 600/610/810
212490	Round punch shelf & bracket set	DKN 900/911/912/912C, DNE 900/910/911, DNF 900/910/911
212094	Round punch shelf & bracket set	DX 300/301/302/302C/312C, DKN 300/301/302/302C/312C
212095	Round punch shelf & bracket set	DX 400/410/401/402/402C/412C, DVS 400/401, DKN 400/401
212808	Round punch shelf & bracket set	DR200
211854	Round punch shelf & bracket set	DG800/850
297071	Round punch shelf	ADP 21/200C/210C
297072	Round punch shelf	ADP 31/300C/310C
212192	Square punch shelf	DP 41/43/43C
212193	Square punch shelf	DP 61/63/63C
211064	Stainless wire shelf & bracket set	DF 611/612 DH 611/612




212246




212266




212490




212094



212095




297071



297072



212192



212193




211064

## Stacking Kit

### Stacking kit and suitable models

Product code	Model	Suitable oven models
212822	OD40	DKN402, DVS402, DX402
212823	OD60	DKN602, DVS602, DX602
212806	ODN26	DNE400/410, DNF400/410 (including models suitable for OD40)
212807	ODN28	DNE600/610, DNF600/610 (including models suitable for OD60)

### OD/ODN



### Important Notes:

- The stacking units must belong to the same series (among 400 series or among 600 series)
- Do not stack 400 series on 600 series. If the upper 400 series is not fixed well, it is easy to topple
- Stacking of old models is forbidden. Lower unit must be new model.

## Earthquake Countermeasure

### Seismic isolation rubber



Product code	296902
Material	Urethane elastomer
Max. load	~100kg (by 4 pcs.)
Size	W50 x D50 x H5mm 1set: 4pcs.


- Pasted at the bottom to prevent unit from falling
- Three layer structure made of urethane elastomer which absorbs 90% impact

## Other Optional Accessories

### Cable port



### N<sub>2</sub> gas inlet device



### External output terminal



### Exhaust duct



Butterfly duct



single exhaust duct



Double duct to suck in more air



Sample oven with single duct



Sample oven with exhaust duct