


Stands


Stands and suitable oven models

Product code	Stand models	Suitable oven models
211856	ON61	DX 400/410/401/402/402C/412C/600/601/602/610/602C/612C, DVS 400/401/402/402C/412C/600/601/602/611/612/602C/612C, DG400 DKM 300C/310C/400C/410C/600C/610C, DKN 300/301/302/302C/312C/400/401/402/402C/412C/600/601/602/611/612/602C/612C, DNE 400/401/410/411/600/610/611, DNF 400/410/600/610
212802	ONS60	DR200
415464	OP43 (without caster)	DF 411/412, DH 411/412
415465	OP63 (without caster)	DF 611/612, DH 611/612
415466	OP46 (with caster, front side with two stoppers)	DF 411/412, DH 411/412
415467	OP66 (with caster, front side with two stoppers)	DF 611/612, DH 611/612


ON61



OP43/63, OP46/66 (casters front side with two stoppers)




ONS60



Shelves

Stainless steel shelf & bracket set and suitable oven models


Product code	Punching shape	Suitable oven models
212246	Round punch shelf & bracket set	DVS 402/402C/412C, DG400, DKN 402/402C/412C, DNE 400/401/410/411, DNF 400/410
212266	Round punch shelf & bracket set	DX 600/601/602/610/602C/612C, DVS 600/601/602/611/612/602C/612C, DKN 600/601/602/611/612/602C/612C/810/811/812/812C, DNE 600/610/611/810/811, DNF 600/610/810
212490	Round punch shelf & bracket set	DKN 900/911/912/912C, DNE 900/910/911, DNF 900/910/911
212094	Round punch shelf & bracket set	DX 300/301/302/302C/312C, DKN 300/301/302/302C/312C
212095	Round punch shelf & bracket set	DX 400/410/401/402/402C/412C, DVS 400/401, DKN 400/401
212808	Round punch shelf & bracket set	DR200
211854	Round punch shelf & bracket set	DG800/850
297071	Round punch shelf	ADP 21/200C/210C
297072	Round punch shelf	ADP 31/300C/310C
212192	Square punch shelf	DP 41/43/43C
212193	Square punch shelf	DP 61/63/63C
211064	Stainless wire shelf & bracket set	DF 611/612 DH 611/612




212246




212266




212490




212094



212095




297071



297072



212192



212193




211064

Stacking Kit

Stacking kit and suitable models

Product code	Model	Suitable oven models
212822	OD40	DKN402, DVS402, DX402
212823	OD60	DKN602, DVS602, DX602
212806	ODN26	DNE400/410, DNF400/410 (including models suitable for OD40)
212807	ODN28	DNE600/610, DNF600/610 (including models suitable for OD60)

OD/ODN



Important Notes:

- The stacking units must belong to the same series (among 400 series or among 600 series)
- Do not stack 400 series on 600 series. If the upper 400 series is not fixed well, it is easy to topple
- Stacking of old models is forbidden. Lower unit must be new model.

Earthquake Countermeasure

Seismic isolation rubber



Product code	296902
Material	Urethane elastomer
Max. load	~100kg (by 4 pcs.)
Size	W50 x D50 x H5mm 1set: 4pcs.


- Pasted at the bottom to prevent unit from falling
- Three layer structure made of urethane elastomer which absorbs 90% impact

Other Optional Accessories

Cable port



N₂ gas inlet device



External output terminal



Exhaust duct



Butterfly duct



single exhaust duct



Double duct to suck in more air



Sample oven with single duct



Sample oven with exhaust duct