

SINCE 1889



SPECIAL OFFER

ROTARY EVAPORATOR BASIC SET

RE-301 ROTARY EVAPORATOR

Yamato's rotary evaporators are built with strong frames making this equipment highly dependable and extremely durable. Yamato rotavaps are known for the steady rotation speeds, high uniformity within the sample, and flexibility of solvent usage.

RE-301AW

LIST PRICE: ~~\$2,860~~

SALE PRICE: \$2,150

GLASSWARE A

standard glass set where condenser is tilted. It is suitable for distillation, concentration and collection of samples.



GLASSWARE B

vertically set, suitable for distillation, concentration, and collection of samples regardless of the boiling point.



RE-301BW

LIST PRICE: ~~\$3,230~~

SALE PRICE: \$2,300



Glassware C (Dry Ice Use)

LIST PRICE: ~~\$1,150~~

SALE PRICE: \$575

50% OFF

Vertically set, suitable for distillation, concentration, and collection of samples which have low boiling points.

SPECIAL PRICE



CF-301 (Optional Item)

LIST PRICE: ~~\$2,570~~

SALE PRICE: \$1,850

- Cooling water circulator
- Closed system with excellent cooling capacity
- Temperature range -20°C to room temp.
- Water capacity ~ 4L (Liquid 3L)
- Refrigerator/Coolant: Air cooling 450W/R404A

RE-301		
Performance	Rotational Speed Range	20 - 250 rpm
Structure	Motor	DC Brushless motor
	Lift	One touch automatic lift
	Material	Diecast Metal
Standard	Power Supply	100 - 240 VAC, 1.5 A
Operational Requirements	Glassware	A, B or C
	Water Bath	BM-500 Water Bath Digital 4.0L 115V 50/60Hz, 12.5A BM-510 Water Bath Digital 4.0L 220V 50/60Hz, 6.5A
CF-301	CF-301 Cooling Water Circulator 4L	115V/220V 50/60Hz, 3.9A/7A

YAMATO SCIENTIFIC AMERICA INC.
925 WALSH AVE. SANTA CLARA, CA 95050

Toll free: 1-800-292-6286
Tel: 1-408-235-7725 / Fax: 1-408-235-7730
customerservice@yamato-usa.com
www.yamato-usa.com

SINCE 1889



SPECIAL OFFER

ROTARY EVAPORATOR BASIC SET

RE-301CW ROTARY EVAPORATOR

Yamato's rotary evaporators are built with strong frames making this equipment highly dependable and extremely durable. Yamato rota-vaps are known for the steady rotation speeds, high uniformity within the sample, and flexibility of solvent usage.

LIST PRICE: ~~\$3,220~~

SALE PRICE: \$2,100



GLASSWARE SET C
vertically set, suitable for distillation, concentration, and collection of samples which have low boiling points.

- Set consists of:
- * RE-301 Rotary Evaporator
 - * 4L water bath
 - * Glassware Set C
 - 1L Evaporating Flask
 - 1L Receiving Flask

RE-301		
Performance	Rotational Speed Range	20 - 250 rpm
Structure	Motor	DC Brushless motor
	Lift	One touch automatic
	Material	Diecast Metal
Standard	Power Supply	100 - 240 VAC, 1.5 A
Operational Requirements	Glassware	C
	Water Bath	BM-500 Water Bath Digital 4.0L 115V 50/60Hz, 12.5A BM-510 Water Bath Digital 4.0L 220V 50/60Hz, 6.5A

YAMATO SCIENTIFIC AMERICA INC.
925 WALSH AVE. SANTA CLARA, CA 95050

Toll free: 1-800-292-6286
Tel: 1-408-235-7725 / Fax: 1-408-235-7730
customerservice@yamato-usa.com
www.yamato-usa.com