


## Stands


### ■ Stands and suitable oven models

Product code	Model	Suitable oven models
211180	ON30	DKM300C/310C, DKN302C/312C, DNF301
211856	ON61	DX400/600 Series, DVS400/600 Series, DG400 Series, DKM300/400/600 Series, DKN300/400/600 Series, DNE400/600 Series, DNF400/600 Series
212348	OT42	DNE401/411, DNF401/411
212349	OT62	DNE601/611, DNF601/611
212477	OH41	DN411IE
212478	OH61	DN611IE
212801	ONS30	DX302C/312C
212802	ONS60	DX402C/412C/602C/612C, DR200
415464	OP43	DF411/412, DH411/412 (stand without caster)
415465	OP63	DF611/612, DH611/612 (stand without caster)
415466	OP46	DF411/412, DH411/412 (stand with caster wheels and stopper in front)
415467	OP66	DF611/612, DH611/612 (stand with caster wheels and stopper in front)


### ■ ON30




### ■ OT42/62




### ■ ON61




### ■ OH41/61



### ■ ONS60




### ■ OP43/63, OP46/66 (casters front side with two stoppers)




## Shelves

### ■ Stainless steel shelf & bracket and suitable oven models

Product code	Description	Model	Suitable oven models
211063	Shelf and bracket set stainless wire (loading up to 30 kg/shelf)	ODQ10	DF/DH412, DN411IE
211064	Shelf and bracket set stainless wire (loading up to 30 kg/shelf)	ODQ20	DF/DH612, DN611IE
211098	Shelf and bracket set stainless punch (loading up to 15kg/shelf)	ODQ30	DF/DH412, DN411IE
211099	Shelf and bracket set stainless punch (loading up to 15 kg/shelf)	ODQ40	DF/DH612, DN611IE
211854	Shelf and bracket set		DG800 Series
212094	Shelf and bracket set stainless punch		300 Series of DX, DKM, DKN, DNF
212095	Shelf and bracket set stainless punch		400 Series of DX, DVS, DKN
212192	Shelf		DP41/43/43C
212193	Shelf		DP61/63/63C
Q110204006	Shelf		DP83C
Q110504007	Shelf		DP103C
212246	Shelf & bracket set stainless punch		400 Series of DVS, DG, DKM, DKN, DNE, DNF
212266	Shelf & bracket set stainless punch		600 Series of DX, DVS, DKM and 600 to 800 Series of DKN, DNE, DNF
212490	Shelf & bracket set stainless punch		900 Series of DKN, DNE, DNF
212678	Stainless steel wire shelf and bracket set. Withstand load up to 30 kg/shelf		DES830/DTS830
212679	Stainless steel punching shelf and bracket set. Withstand load up to 15 kg/shelf	ODE50	DES830/DTS830
212808	Shelf		DR200
212919	Basket-type shelf placed on top of stD shelves. Withstand load up to 15kg/shelf)	ODE12	DES830/DTS830
212924	Shelf and bracket set stainless mesh (loading up to 15 kg/shelf, 30mm deep/ designed to be stacked on std stainless wire shelves)	ODT12	DF/DH412
212925	Shelf and bracket set stainless mesh (loading up to 15 kg/shelf, 30 mm deep/ designed to be stacked on standard stainless wire shelves)	ODT14	DF/DH612
297071	Shelf		ADP21/200C/210C
297072	Shelf		ADP31/300C/310C
YSA0000071	Shelf and bracket set		DF/DH832




212246




212266




212490




212094



212095



297071



297072



212192



212193



211064

## Stacking Kit


### ■ Stacking kit and suitable models

Product code	Model	Suitable oven models
212803	ODK80	DX302C/312C
212804	ODK82	DX402C/412C
212805	ODK84	DX602C/612C
212806	ODN26	400 Series of DNE, DNF (including models suitable for OD40)
212807	ODN28	600 Series of DNE, DNF (including models suitable for OD60)
212822	OD40	400 Series of DX, DVS, DKM, DKN
212823	OD60	600 Series of DX, DVS, DKM, DKN
213700	ODF48	DF412/612, DH412/612
281458	ODM44	DNF301



### ■ Important Notes:

- The stacking units must belong to the same series (among 400 series or among 600 series)
- Do not stack 400 series on 600 series. If the upper 400 series is not fixed well, it is easy to topple
- Stacking of old models is forbidden. Lower unit must be new model.

### ■ OD/ODN



### ■ ODK



## Earthquake Countermeasure

### Seismic isolation rubber




Product code	296902
Material	Urethane elastomer
Max. load	~100kg (by 4 pcs.)
Size	W50 x D50 x H5mm 1set: 4pcs.

- Pasted at the bottom to prevent unit from falling
- Three layer structure made of urethane elastomer which absorbs 90% impact

## Other Optional Accessories

### Cable port



### Exhaust duct



Butterfly duct



Double duct to suck in more aire



Exhaust duct for DF and DH 213703 / 213704

### N<sub>2</sub> gas inlet device



### External output terminal

