


Stands


Stands and suitable oven models

| Product code | Model | Suitable oven models |
|--------------|-------|---|
| 211180 | ON30 | DKM300C/310C, DKN302C/312C, DNF301 |
| 211856 | ON61 | DX400/600 Series, DVS400/600 Series, DG400 Series, DKM300/400/600 Series, DKN300/400/600 Series, DNE400/600 Series, DNF400/600 Series |
| 212348 | OT42 | DNE401/411, DNF401/411 |
| 212349 | OT62 | DNE601/611, DNF601/611 |
| 212477 | OH41 | DN411IE |
| 212478 | OH61 | DN611IE |
| 212801 | ONS30 | DX302C/312C |
| 212802 | ONS60 | DX402C/412C/602C/612C, DR200 |
| 415464 | OP43 | DF411/412, DH411/412 (stand without caster) |
| 415465 | OP63 | DF611/612, DH611/612 (stand without caster) |
| 415466 | OP46 | DF411/412, DH411/412 (stand with caster wheels and stopper in front) |
| 415467 | OP66 | DF611/612, DH611/612 (stand with caster wheels and stopper in front) |


ON30




OT42/62




ON61




OH41/61



ONS60




OP43/63, OP46/66 (casters front side with two stoppers)




Shelves

■ Stainless steel shelf & bracket and suitable oven models


| Product code | Description | Model | Suitable oven models |
|--------------|---|-------|---|
| 211063 | Shelf and bracket set stainless wire (loading up to 30 kg/shelf) | ODQ10 | DF/DH412, DN411IE |
| 211064 | Shelf and bracket set stainless wire (loading up to 30 kg/shelf) | ODQ20 | DF/DH612, DN611IE |
| 211098 | Shelf and bracket set stainless punch (loading up to 15kg/shelf) | ODQ30 | DF/DH412, DN411IE |
| 211099 | Shelf and bracket set stainless punch (loading up to 15 kg/shelf) | ODQ40 | DF/DH612, DN611IE |
| 211854 | Shelf and bracket set | | DG800 Series |
| 212094 | Shelf and bracket set stainless punch | | 300 Series of DX, DKM, DKN, DNF |
| 212095 | Shelf and bracket set stainless punch | | 400 Series of DX, DVS, DKN |
| 212192 | Shelf | | DP41/43/43C |
| 212193 | Shelf | | DP61/63/63C |
| Q110204006 | Shelf | | DP83C |
| Q110504007 | Shelf | | DP103C |
| 212246 | Shelf & bracket set stainless punch | | 400 Series of DVS, DG, DKM, DKN, DNE, DNF |
| 212266 | Shelf & bracket set stainless punch | | 600 Series of DX, DVS, DKM and 600 to 800 Series of DKN, DNE, DNF |
| 212490 | Shelf & bracket set stainless punch | | 900 Series of DKN, DNE, DNF |
| 212678 | Stainless steel wire shelf and bracket set. Withstand load up to 30 kg/shelf | | DES830/DTS830 |
| 212679 | Stainless steel punching shelf and bracket set. Withstand load up to 15 kg/shelf | ODE50 | DES830/DTS830 |
| 212808 | Shelf | | DR200 |
| 212919 | Basket-type shelf placed on top of stD shelves. Withstand load up to 15kg/shelf) | ODE12 | DES830/DTS830 |
| 212924 | Shelf and bracket set stainless mesh (loading up to 15 kg/shelf, 30mm deep/ designed to be stacked on std stainless wire shelves) | ODT12 | DF/DH412 |
| 212925 | Shelf and bracket set stainless mesh (loading up to 15 kg/shelf, 30 mm deep/ designed to be stacked on standard stainless wire shelves) | ODT14 | DF/DH612 |
| 297071 | Shelf | | ADP21/200C/210C |
| 297072 | Shelf | | ADP31/300C/310C |
| YSA0000071 | Shelf and bracket set | | DF/DH832 |




212246




212266




212490




212094



212095




297071



297072



212192



212193



211064

Stacking Kit


■ Stacking kit and suitable models

| Product code | Model | Suitable oven models |
|--------------|-------|---|
| 212803 | ODK80 | DX302C/312C |
| 212804 | ODK82 | DX402C/412C |
| 212805 | ODK84 | DX602C/612C |
| 212806 | ODN26 | 400 Series of DNE, DNF (including models suitable for OD40) |
| 212807 | ODN28 | 600 Series of DNE, DNF (including models suitable for OD60) |
| 212822 | OD40 | 400 Series of DX, DVS, DKM, DKN |
| 212823 | OD60 | 600 Series of DX, DVS, DKM, DKN |
| 213700 | ODF48 | DF412/612, DH412/612 |
| 281458 | ODM44 | DNF301 |



■ Important Notes:

- The stacking units must belong to the same series (among 400 series or among 600 series)
- Do not stack 400 series on 600 series. If the upper 400 series is not fixed well, it is easy to topple
- Stacking of old models is forbidden. Lower unit must be new model.

■ OD/ODN



■ ODK



Earthquake Countermeasure

Seismic isolation rubber




| | |
|--------------|---------------------------------|
| Product code | 296902 |
| Material | Urethane elastomer |
| Max. load | ~100kg (by 4 pcs.) |
| Size | W50 x D50 x H5mm 1set: 4pcs. |

- Pasted at the bottom to prevent unit from falling
- Three layer structure made of urethane elastomer which absorbs 90% impact

Other Optional Accessories

Cable port



Exhaust duct



Butterfly duct



Double duct to suck in more aire



Exhaust duct for DF and DH 213703 / 213704

N₂ gas inlet device



External output terminal

