

SINCE 1889



# BENCHTOP VACUUM DRYING OVEN

## ADP200C/210C/300C/310C

**Standard small size benchtop vacuum drying oven with enhanced safety features**

### Features

- Easy input of parameters and settings
- Supports Auto start and Auto quick stop
- Silicon rubber door seal prevents air from leaking
- Self diagnostic and overheating prevention functions
- Independent overheating prevention device for each circuit
- Customizable with N<sub>2</sub> gas inlet and communication ports
- Easy maintenance



Control Panel



### Specifications

Model	ADP200C	ADP210C	ADP300C	ADP310C
Operating temp. range	40~240°C			
Operating pressure range	101 to 0.1 kPa (760 to 1 Torr) at absolute pressure			
Temp. control accuracy	±1.5°C (at 240 °C)			
Max. temp. reaching time	~70 min.		~100 min.	
Heater/Heating method	Mica heater, 0.68kW/ Decompressed chamber direct heating		Mica heater, 1.05kW / Decompressed chamber direct heating	
Internal dimensions(WxDxH mm)	200×250×200		300×300×300	
Internal capacity	10L		27L	
Shelves	2pcs		3pcs	
Power source 50/60Hz	AC115V 6A	AC220V 3.5A	AC115V 9.5A	AC220V 5A
Weight	~30kg		~55kg	