

ROTARY EVAPORATOR OVERVIEW



RE202-A / 212-A

REV202M-A / 212M-A

- Market leader in rotation speed range at 5~315 rpm
- Designed with set inversion function
- Standard equipped with condenser insulation kit
- Convenient height adjustment
- Integration of servo type motor with strong high torque
- Cables stored neatly, simply and smartly
- Option for 2 types of bath: 5L water or oil bath 115~240VAC
- Option for a vacuum controller



RE301

- Easy manual operation
- Digitally set and displayed motor speed (rpm)
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC



RE601

- Programmable
- Digital setting and display for motor speed (rpm), vacuum control and vapor temperature
- Equipped with vacuum regulator and vapor temperature indicator
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC



RE801

- Programmable
- Bath synchronized control
- Automatic distillation and 53 solvent presets
- Digital setting and display for motor speed (rpm), vacuum control and vapor temperature
- Equipped with vacuum regulator and vapor temperature indicator
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC